



CENTRO DE
VIGILANCIA SANITARIA VETERINARIA
Universidad Complutense



PEN SIDE TEST FOR AFRICAN SWINE FEVER ANTIBODY DETECTION

jmvizcaino@visavet.ucm.es
Av/ Puerta de Hierro s/n.
28040 Madrid.

Tel: (34) 913944082
Fax: (34) 913943908

PEN SIDE TEST FOR AFRICAN SWINE FEVER ANTIBODY DETECTION

1. MATERIALS

- Single channel pipettes 10-100 μ l
- Latex or nitrile gloves
- Pen side test : 11.PPA.K.4: INGEZIM PPA CROM Antibody detection (INGENASA [®])

2. METHODOLOGY

2.1) Sample dilution: The serum sample needs to be at 1/50 dilution. Take 10 μ l of the serum sample and put it into the diluent tube. Mix them properly.

2. 2) Sample addition: add 150 μ l of the diluted sample in the round window of the device.

2.3) Wait 10 minutes for reading the results. No result must be considered after 10 minutes


3. INTERPRETATION OF THE RESULTS

Test validation: a blue line must always appear in control area (C) to ensure the correct working of the test. If there is no blue line in control area (C) the assay cannot be considered valid.

Result interpretation:

- POSITIVE SAMPLE: red/pink line in the test area (T)
- NEGATIVE SAMPLE: no red/pink line in the test area (T)

*Positive samples should be confirmed by techniques recommended by OIE





For more information, please contact us:

jmvizcaino@visavet.ucm.es

estefania.cadenas@ucm.es

VISAVET Center and Department of Animal Health
Facultad de Veterinaria. Universidad Complutense de Madrid
Av/ Puerta de Hierro s/n.
28040 Madrid, Spain

Tel: (34) 913944082

Fax: (34) 913943908